




## VITREX<sup>®</sup> Safety Containers

- ◆ VITREX<sup>®</sup> Safety Containers are designed for hospitals and health care institutions for clinical disposal such as lancets, needles, scalpels, syringes and other types of sharps.
- ◆ USER SAFETY - designed to meet UN requirement no. 3291.
- ◆ USER CONVENIENCE - Each type of lid is specially designed to help the user dispose of needles, syringes, etc.
- ◆ Lid with reversible closure can be safely fixed.
- ◆ The Safety Containers are available in all common sizes.
- ◆ Environmentally friendly (polypropylene) materials are used.

## PRODUCTS

SIZE	DIAMETER TOP / HEIGHT	ARTICLE NO.	PIECES PER CARTON	PIECES PER PALLET
0.1 L	4.7 / 7.5 cm.	193700	100 pcs.	7,400 pcs.
0.2 L	4.7 / 12.9 cm.	193800	100 pcs.	6,400 pcs.
0.5 L	9.0 / 10.8 cm.	194000	100 pcs.	3,200 pcs.
1.5 L	13.7 / 12.5 cm.	194200	100 pcs.	2,000 pcs.
2.0 L	13.7 / 18.0 cm.	194300	100 pcs.	1,600 pcs.
5.0 L	21.1 / 19.5 cm.	194700	100 pcs.	405 pcs.
11.0 L	26.0 / 25.5 cm.	195000	100 pcs.	198 pcs.
21.0 L	29.0 / 39.5 cm.	195300	100 pcs.	180 pcs.

Syringes



Insulin pens  
(different brands/  
diameters)

Scalpels and  
other types of  
sharps

Designed for disposal of syringes, insulin pen-needles, lancets, scalpels and other types of sharps



Lid for 0.5 L to 21,0 L

Shows various functions for disposal.  
(Here used for disposal of syringe).



Lid for 0,1 L and 0,2 L.

Designed especially for Home Care nursing and doctors.

All Safety Containers meet UN requirement no 3291 which a.o. means that you can drop the safety container from 1,5 meter without breakage or leakage