



VITREX[®] Soft BLOOD GLUCOSE LANCETS 23, 28, 30 and 32G

- ◆ RANGE OF BLOOD VOLUMES 23, 28, 30 and 32G lancets
- ◆ GOOD BLOOD FLOW triple ground needle point
- ◆ VERSATILE fit most commonly used lancing devices
- ◆ STERILE sterilised by gamma irradiation at minimum 25kGy
- ◆ CONVENIENT boxes have dispensing feature

PRODUCTS

DESCRIPTION	GAUGE SIZE	PIECES/BOX	ARTICLE NUMBER
Sterile blood lancet	23G (0.65mm)	200	3512-13
Sterile blood lancet	28G(0.38mm)	200	3522-13
Sterile blood lancet	30G (0.31mm)	200	3542-13
Sterile blood lancet	32G (0.26mm)	200	3572-13

PACKAGING

PRIMARY BOX	CASE	SHIPPER	PALLET
200 pcs	10 x 200 2,000 pcs	10 x cases 20,000 pcs	480.000 pcs

PRODUCT SPECIFICATION

Needle

MATERIAL	Stainless Steel Grade 304/1 Cr18Ni9		
EXPOSED LENGTH	3.3mm (+/- 0.33mm)		
POINT	3 facet (hypodermic) profile		
NOMINAL DIAMETER	0.65mm (23G)		
	0.38mm (28G)		
	0.31mm (30G)		
	0.26mm (32G)		

Moulded Body/Cap

MATERIAL	LDPE (LOW DENSITY POLYETHYLENE)		
OVERALL LENGTH	32.5mm		
COLOUR	Light blue		23G
	Red		28G
	Green		30G
	Dark blue		32G

Sterility

METHOD	Gamma irradiation
MINIMUM DOSAGE	25kGY
SHELF LIFE	5 years

VITREX® SoFT IS DESIGNED FOR USE WITH MOST POPULAR TYPES OF FINGER LANCING DEVICE