



VITREX[®] Steel GENERAL PURPOSE CAPILLARY BLOOD LANCET

- ◆ LOW COST simple, flat, steel lancet
- ◆ SAMPLE VOLUME excellent blood flow
- ◆ HIGH QUALITY hardened surgical steel
- ◆ MINIMUM DISCOMFORT modern asymetric design
- ◆ STERILE sterilised by gamma radiation to DIN EN 552/556

PRODUCT

PRODUCT	DESCRIPTION	ARTICLE NUMBER
VITREX® Steel	Flat blade sterile lancet	3502-13

PACKAGING

PRIMARY BOX	CASE	SHIPPER	PALLET
200 pcs	50 x Primaries (10,000pcs)	4 x cases (40,000pcs)	720.000 pcs

PRODUCT SPECIFICATION

Blade

MATERIAL	Stainless Steel Grade X 5 Cr18Ni9 to DIN 17440
TOTAL LENGTH	49.00mm (+/- 0.50mm)
EXPOSED LENGTH	3.0mm (+/- 0.30mm)
POINT	12° asymmetric point with 18° angle

Sterility

METHOD	Gamma irradiation by cobalt 60 source to DIN EN 552/556
SHELF LIFE	5 years
WRAPPING PAPER	
TOP (PRINTED)	Supraferm 40g/m ² kraft paper 20g/m ² polyethylene coated
BOTTOM (TRANSPARENT)	Supraferm 30g/m ² bleached pergamyne 20g/m ² polyethylene coated