



VITREX® STERILANCE® Flex STERILE SAFETY LANCET

- ◆ USER SAFETY is assured because the needle is fully shielded before and after use
- ◆ COMFORT is provided by high quality tri-bevelled needle with high-speed penetration
- ◆ USER CONVENIENCE - adjustable design provides three different penetration depths from a single unit - easy to use pre-loaded and pressure activated device
- ◆ PATIENT SAFETY sterile, single-use, unable to be re-used
- ◆ PRESSURE ACTIVATION ensures a consistent blood flow with every use
- ◆ APPLICATION - can be used for blood glucose and other diagnostic tests

COLOUR	GAUGE/DEPTH	PIECES/PACK	ARTICLE NO	DESCRIPTION
Blue	26G / 1.2 mm	100 pcs	399913	Adjustable Safety Lancet
	26G / 1.8 mm			
	26G / 2.4 mm			

Shelf life 3 years



- ▶ Carefully pull off the protective cap



- ▶ **WARNING:** Do not use if protective cap has previously been removed

- ▶ Adjust the adjustable tip to choose the suitable depth



- ▶ Push the Adjustable Safety Lancet firmly onto the chosen site to fire



- ▶ Dispose of the Adjustable Safety Lancet in a suitable safety container

PACKAGING

Primary pack:	100 pcs
Shipper:	10x100 primaries
Pieces per shipper:	1000 pcs
Pallet	84.000 pcs