




VITREX[®] Safety Containers

- ◆ VITREX[®] Safety Containers are designed for hospitals and health care institutions for clinical disposal such as lancets, needles, scalpels, syringes and other types of sharps.
- ◆ USER SAFETY - designed to meet UN requirement no. 3291.
- ◆ USER CONVENIENCE - Each type of lid is specially designed to help the user dispose of needles, syringes, etc.
- ◆ Lid with reversible closure can be safely fixed.
- ◆ The Safety Containers are available in all common sizes.
- ◆ Environmentally friendly (polypropylene) materials are used.

PRODUCTS

SIZE	DIAMETER TOP/HEIGHT	ARTICLE NO.	PCS PER CARTON	PCS PER PALLET
0,1 L	4,7 / 7,5 CM	193700	100 PCS	8400 PCS
0,2 L	4,7 / 12,9 CM	193800	100 PCS	6400 PCS
0,5 L	9,0 / 10,8 CM	194000	100 PCS	3200 PCS
1,5 L	13,7 / 12,5 CM	194200	100 PCS	2000 PCS
2,0 L	13,7 / 18,0 CM	194300	100 PCS	1600 PCS
5,0 L	21,1 / 19,5 CM	194700	N/A	450 PCS
11,0 L	26,0 / 25,5 CM	195000	N/A	198 PCS
21,0 L	29,0 / 39,5 CM	195300	N/A	180 PCS

Syringes



Insulin pens
(different brands/

Scalpels and
other types of
sharps

- ▶ Designed for disposal of syringes, insulin pen-needles, lancets, scalpels and other types of sharps



- ▶ Lid for 0.5 L to 21,0 L

Shows various functions for disposal.
(Here used for disposal of syringe).



- ▶ Lid for 0,1 L and 0,2 L.

Designed especially for Home Care nursing and doctors.

All Safety Containers meet UN requirement no 3291 which a.o. means that you can drop the safety container from 1,5 meter without breakage or leakage