



DIABETES

## Chirana® Insulin and Tuberculin Syringes

Chirana Insulin and Tuberculin syringes are with integrated needle for single use. The needle cap is made of LATEX-FREE rubber gasket. The syringes are individually packed in boxes of 100 pcs. They are pyrogen-free and non-toxic.

**USER CONVENIENCE:** Needle permanently fixed to syringe body.

**USER SAFETY:** Integrated needle cap.

**COMFORT:** Grinded, precision lancet point and silicone treated tip for smooth penetration.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO. Dead space is minimized to reduce waste and transmission of infectious agents. The dead space does not exceed the limits (EN ISO 8537:2016, art. 5.11.1).

**FDA CLEARED \***



**Packaging:**  
Primary packaging: 100 pcs.  
Secondary packaging: 2,000 pcs.  
Pallet: 48,000 pcs.

Transit packaging

Pallet

\* FDA clearance for Insulin syringes.

### Technical Specifications

| Type          | Vol. ml. | Needle      | Size Metric | Article No: |
|---------------|----------|-------------|-------------|-------------|
| Tuberculin    | 1        | 29G x 1/2"  | 12.7 mm     | CHTUB0129   |
| Tuberculin    | 1        | Luer        |             | CHTUB01     |
| Insulin U 40  | 1        | 29G x 1/2"  | 12.7 mm     | CHINS4129   |
| Insulin U 100 | 1        | 27G x 1/2"  | 12.7 mm     | CHINS0127   |
| Insulin U 100 | 1        | 29G x 5/16" | 8 mm        | CHINS01298  |
| Insulin U 100 | 1        | 29G x 1/2"  | 12.7 mm     | CHINS0129   |
| Insulin U 100 | 1        | 30G x 1/2"  | 12.7 mm     | CHINS0130   |
| Insulin U 100 | 1        | 30G x 5/16" | 8 mm        | CHINS01308  |
| Insulin U 100 | 1        | 31G x 4/16" | 6 mm        | CHINS01316  |
| Insulin U 100 | 1        | 31G x 5/16" | 8 mm        | CHINS01318  |
| Insulin U 100 | 1        | Luer        |             | CHINS01     |
| Insulin U 100 | 0.5      | 29G x 1/2"  | 12.7 mm     | CHINS00529  |
| Insulin U 100 | 0.5      | 30G x 1/2"  | 12.7 mm     | CHINS00530  |
| Insulin U 100 | 0.5      | 30G x 5/16" | 8 mm        | CHINS005308 |
| Insulin U 100 | 0.5      | 31G x 5/16" | 8 mm        | CHINS005318 |
| Insulin U 100 | 0.5      | 31G x 4/16" | 6 mm        | CHINS005316 |
| Insulin U 100 | 0.3      | 29G x 1/2"  | 12.7 mm     | CHINS00329  |
| Insulin U 100 | 0.3      | 30G x 1/2"  | 12.7 mm     | CHINS00330  |
| Insulin U 100 | 0.3      | 31G x 5/16" | 8 mm        | CHINS003318 |
| Insulin U 100 | 0.3      | 31G x 4/16" | 6 mm        | CHINS003316 |

CHIRANA Syringes are available in many sizes.  
Ask your diabetes educator or pharmacist for advice.

**Shelf life 5 years**