



## Chirana® 3-part disposable syringes and Chirana® 3-part disposable syringes with needles

Chirana 3-part disposable syringes without and with needle are made of polypropylene, polyethylene and lubricated synthetic rubber. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The syringes are available with Luer and Luer lock. The syringe with needle is available in other needle combinations.

**USER CONVENIENCE:** Ergonomic design, smooth flow, highly transparent barrel, tight, black extended graduation and easy to read.

**USER SAFETY:** Individually blister packed and the needle is side-packed in a single blister package.

**COMFORT:** Designed tip for best connection, scratch-free for glove manipulation and centrally positioned Luer lock tip.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO. Low dead space volume for 1 ml.

### Chirana 3-part disposable syringes

Vol. (ml)	Pcs./Box	Pcs./Carton	Article No. LUER	Article No. LUER LOCK
1	100	2000	CH03001L	CH03001LL
2	100	3000	CH03002L	CH03002LL
3	100	3000	CH03003L	CH03003LL
5	100	2400	CH03005L	CH03005LL
10	100	1600	CH03010L	CH03010LL
20	50	800	CH03020L	CH03020LL
30	50	400	CH03030L	CH03030LL
50		100	CH03050L	CH03050LL

Shelf life 5 years

### Chirana 3-part disposable syringes with needles

Specification	Pcs./Box	Pcs./Carton	Article No.
2 ml with needle 0.60 x 30 mm (23G x 1 1/4")	100	3000	CH03002L0630
5 ml with needle 0.70 x 30 mm (22G x 1 1/4")	100	2400	CH03005L0730
10 ml with needle 0.80 x 40 mm (21G x 1 1/2")	100	1600	CH03010L0840
20 ml with needle 0.80 x 40 mm (21G x 1 1/2")	50	800	CH03020L0840
50 ml Luer with needle 1.20 x 40 mm (18G x 1 1/2")	50	450	CH03050L1240

Shelf life 5 years