



Medoject® Hypodermic needles

Colour Code		Size Metric	Size G/°	Article No.
Light yellow	●	0.30 x 12 mm	30G x 1/2"	CH30012
Grey	●	0.40 x 12 mm	27G x 1/2"	CH27012
Grey	●	0.40 x 20 mm	27G x 4/5"	CH27045
Grey	●	0.40 x 20 mm	27G x 3/4"	CH27034
Grey	●	0.40 x 40 mm	27G x 1 1/2"	CH27112
Brown	●	0.45 x 12 mm	26G x 1/2"	CH26012
Brown	●	0.45 x 16 mm	26G x 5/8"	CH26058
Brown	●	0.45 x 25 mm	26G x 1"	CH26100
Orange	●	0.50 x 16 mm	25G x 5/8"	CH25058
Orange	●	0.50 x 20 mm	25G x 4/5"	CH25045
Orange	●	0.50 x 25 mm	25G x 1"	CH25100
Orange	●	0.50 x 40 mm	25G x 1 1/2"	CH25112
Purple	●	0.55 x 25 mm	24G x 1"	CH24100
Purple	●	0.55 x 38 mm	24G x 1 1/2"	CH24112
Blue	●	0.60 x 16 mm	23G x 5/8"	CH23058
Blue	●	0.60 x 25 mm	23G x 1"	CH23100
Blue	●	0.60 x 30 mm	23G x 1 1/4"	CH23114
Blue	●	0.60 x 40 mm	23G x 1 1/2"	CH23112
Blue	●	0.60 x 60 mm	23G x 2 1/4"	CH23214
Black	●	0.70 x 25 mm	22G x 1"	CH22100
Black	●	0.70 x 30 mm	22G x 1 1/4"	CH22114
Black	●	0.70 x 35 mm	22G x 1 2/5"	CH22125
Black	●	0.70 x 40 mm	22G x 1 1/2"	CH22112
Black	●	0.70 x 50 mm	22G x 2"	CH22200
Green	●	0.80 x 16 mm	21G x 5/8"	CH21058
Green	●	0.80 x 25 mm	21G x 1"	CH21100
Green	●	0.80 x 30 mm	21G x 1 1/4"	CH21114
Green	●	0.80 x 40 mm	21G x 1 1/2"	CH21112
Green	●	0.80 x 50 mm	21G x 2"	CH21200
Yellow	●	0.90 x 25 mm	20G x 1"	CH20100
Yellow	●	0.90 x 30 mm	20G x 1 1/4"	CH20114
Yellow	●	0.90 x 40 mm	20G x 1 1/2"	CH20112
Yellow	●	0.90 x 50 mm	20G x 2"	CH20200
Ivory	○	1.10 x 25 mm	19G x 1"	CH19100
Ivory	○	1.10 x 30 mm	19G x 1 1/4"	CH19114
Ivory	○	1.10 x 40 mm	19G x 1 1/2"	CH19112
Ivory	○	1.10 x 50 mm	19G x 2"	CH19200
Pink	●	1.20 x 25 mm	18G x 1"	CH18100
Pink	●	1.20 x 40 mm	18G x 1 1/2"	CH18112
Pink	●	1.20 x 50 mm	18G x 2"	CH18200
Red purple	●	1.50 x 50 mm	17G x 2"	CH17200
White	○	1.60 x 25 mm	16G x 1"	CH16100
White	○	1.60 x 40 mm	16G x 1 1/2"	CH16112
Blue grey	●	1.80 x 40 mm	15G x 1 1/2"	CH15112

Medoject®

Medoject needles are made of polypropylene hub and stainless steel tube. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The needles are intended for intramuscular, subcutaneous, intravenous and intra-arterial application. They are individually packed in boxes of 100 pcs, cartons of 4000 pcs.

USER CONVENIENCE: Colour-coded transparent hub for flow visual control, metric and gauge measurements on packaging and mainly thin wall tube material.

USER SAFETY: Safety needle is available.

COMFORT: Grinded, precision lancet point and silicone treated for smooth penetration.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Medoject® Hypodermic needles

Medoject Hypodermic needles are intended for intramuscular, subcutaneous, intravenous and intra-arterial application. There is possibility of short bevel lancet point grinding.

Shelf life 5 years