

CHIRAFLON monofilament

CHIRAFLON monofilament is a non-absorbable suture made of PVDF (similar to teflon). Physiologically this thread is completely inert.

Chiraflon monofilament offers optimal smoothness and superior knotting properties. The thread is characterised by its high tensile strength. No splitting occurs even under considerable strain.



Indication: All surgical specialities, especially Orthopaedic, Traumatology (tendons), Plastic and Ophthalmic surgery.

Contraindication: Chiraflon monofilament must not penetrate into the inner parts of the digestive and urinary tracts.

Chemical Composition: Polyvinylidene Fluoride (PVDF)

Sizes: 8/0 to 3 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue

CHIRAFLON monofilament

UK04
136



Non-absorbable - Chiraflon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 45 cm	Blue	TF 6202	
HR 18 	1	5/0	1 x 45 cm	Blue	TF 6211	
	1.5	4/0	1 x 45 cm	Blue	TF 6214	
	1	5/0	1 x 75 cm	Blue	TF 6212	
	1.5	4/0	1 x 75 cm	Blue	TF 6213	
2 x HR 18 	1	5/0	1 x 75 cm	Blue	TF 6216	
	1.5	4/0	1 x 75 cm	Blue	TF 6217	
HR 27 	3	2/0	1 x 75 cm	Blue	TF 6221	
HRX: ½ circle, round body needle - extra strong						
HRX 40 	3.5	0	1 x 90 cm	Blue	TF 6251	
HRX 50 	4	1	1 x 150 cm	Blue	TF 4540	Loop
	3.5	0	1 x 200 cm	Blue	TF 5035	Loop
	4	1	1 x 200 cm	Blue	TF 5040	Loop
HRT: ½ circle, round body needle - tapercut						
HRT 13 	0.7	6/0	1 x 45 cm	Blue	TF 6421	
2 x HRT 17 	1	5/0	1 x 75 cm	Blue	TF 6332	
DRT: 3/8 circle, round body needle - tapercut						
2 x DRT 10 	0.7	6/0	1 x 75 cm	Blue	TF 6301	
2 x DRT 13 	0.7	6/0	1 x 75 cm	Blue	TF 6311	
	1	5/0	1 x 75 cm	Blue	TF 6312	

CHIRAFLON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
GR: Straight, round body needle						
GR 20 	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Blue Blue	TF 6411 TF 6412	
GR 60 	2	3/0	1 x 75 cm	Blue	TF 6462	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Blue Blue	TF 6101 TF 6102	
DS 15 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Blue Blue	TF 6105 TF 6106	
DS 19 	1	5/0	1 x 45 cm	Blue	TF 6110	
	1.5	4/0	1 x 45 cm	Blue	TF 6115	
	2	3/0	1 x 45 cm	Blue	TF 6116	
	3	2/0	1 x 45 cm	Blue	TF 6117	
	1.5	4/0	1 x 75 cm	Blue	TF 6111	
	2	3/0	1 x 75 cm	Blue	TF 6112	
	3	2/0	1 x 75 cm	Blue	TF 6113	
DS 25 	2	3/0	1 x 45 cm	Blue	TF 6126	
	3.5	0	1 x 45 cm	Blue	TF 6128	
	2	3/0	1 x 75 cm	Blue	TF 6122	
	3	2/0	1 x 75 cm	Blue	TF 6123	
DS 30 	3	2/0	1 x 75 cm	Blue	TF 6131	
	3.5	0	1 x 75 cm	Blue	TF 6132	
	4	1	1 x 75 cm	Blue	TF 6133	
DS 38 	2	3/0	1 x 75 cm	Blue	TF 6140	



CHIRAFLON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7 1 1.5	6/0 5/0 4/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Blue Blue Blue	TF 6001 TF 6002 TF 6003	Advanced cutting needle
DSM 15 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Blue Blue	TF 6005 TF 6006	Advanced cutting needle
	1	5/0	1 x 75 cm	Blue	TF 6008	
DSM 19 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Blue Blue	TF 6011 TF 6012	Advanced cutting needle
GS: Straight, reverse cutting needle						
GS 60 	3	2/0	1 x 75 cm	Blue	TF 6167	