





Polydox
monofilament



POLYDOX monofilament

POLYDOX monofilament is a long-term absorbable, synthetic suture made of polydioxanone. The suture offers good handling combined with great strength due to its slow degradation profile.

Polydiox monofilament is ideal for all indications, where an extended wound support of more than 4 weeks is desirable. The threads are available in violet and natural.

The suture is very flexible, pliable and easy to knot. In tissue, the thread can cause a slight initial inflammatory reaction.



Indication: General soft tissue approximation, Orthopaedics (e.g. Joint capsule, Tendons), Fascia Closure, Urology, Gynaecology and Plastic surgery.

Contraindication: Polydiox monofilament should not be used where prolonged approximation of tissues under stress is required.

Sizes: 7/0 to 4 USP

Structure: Monofilament

Absorption method: Hydrolysis

Strength Retention Profile: LONG-TERM HEALING (90 - 98 days)

1 Day	100%
45 Days	50%

Mass Absorption profile: 180 - 240 days

POLYDOX monofilament

UK07
122

Absorbable - Polydiox monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	0.7 1	6/0 5/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0246 DX 0247	
HR 13 	1	5/0	1 x 45 cm	Violet	DX 0248	
	0.7 1 1.5	6/0 5/0 4/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0249 DX 0250 DX 0251	
2 x HR 13 	1	5/0	1 x 75 cm	Violet	DX 5245	
HR 17 	1	5/0	1 x 75 cm	Violet	DX 0254	
HR 18 	1 1.5 2	5/0 4/0 3/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0257 DX 0258 DX 0259	
2 x HR 18 	1	5/0	1 x 75 cm	Violet	DX 5250	
	1.5	4/0	1 x 90 cm	Violet	DX 5251	
HR 20 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0360 DX 0256	
HR 22 	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0261 DX 0262 DX 0263	
2 x HR 22 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 5261 DX 5262	
	2	3/0	1 x 90 cm	Violet	DX 5264	
HR 26 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0366 DX 0365	
2 x HR 26 	3	2/0	1 x 75 cm	Violet	DX 5269	
	1.5 2	4/0 3/0	1 x 90 cm 1 x 90 cm	Violet Violet	DX 5271 DX 5270	

POLYDOX monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 27	1.5	4/0	1 x 75 cm	Violet	DX 0265	
	2	3/0	1 x 75 cm	Violet	DX 0266	
	3	2/0	1 x 75 cm	Violet	DX 0267	
	3.5	0	1 x 75 cm	Violet	DX 0268	
	4	1	1 x 75 cm	Violet	DX 0369	
 2 x HR 27	1.5	4/0	1 x 75 cm	Violet	DX 5265	
	2	3/0	1 x 75 cm	Violet	DX 5266	
	3	2/0	1 x 75 cm	Violet	DX 5267	
	1.5	4/0	1 x 90 cm	Violet	DX 5295	
	2	3/0	1 x 120 cm	Violet	DX 5297	
 HR 30	1.5	4/0	1 x 75 cm	Violet	DX 0269	
	2	3/0	1 x 75 cm	Violet	DX 0270	
	3	2/0	1 x 75 cm	Violet	DX 0271	
	3.5	0	1 x 75 cm	Violet	DX 0272	
	4	1	1 x 75 cm	Violet	DX 0273	
 2 x HR 30	2	3/0	1 x 90 cm	Violet	DX 5275	
 HR 37	2	3/0	1 x 75 cm	Violet	DX 0264	
	3.5	0	1 x 75 cm	Violet	DX 0356	
 HR 40	3	2/0	1 x 75 cm	Violet	DX 0275	
	3.5	0	1 x 75 cm	Violet	DX 0276	
	4	1	1 x 75 cm	Violet	DX 0277	
	3	2/0	1 x 90 cm	Violet	DX 0375	
	4	1	1 x 90 cm	Violet	DX 0377	
 HR 45	4	1	1 x 150 cm	Violet	DX 0279	
	5	2	1 x 90 cm	Violet	DX 0378	



POLYDOX monofilament

UK07
124

Absorbable - Polydox monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 50 	4	1	1 x 90 cm	Violet	DX 0380	
HR 65 	4	1	1 x 150 cm	Violet	DX 4547	Loop
	4	1	1 x 250 cm	Violet	DX 4546	Loop
HR s: ½ circle, round body needle - slim						
HR s 18 	1.5	4/0	1 x 75 cm	Violet	DX 1256	Slim needle
HR s 22 	2	3/0	1 x 75 cm	Violet	DX 1262	Slim needle
HR s 27 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 1272 DX 1273	Slim needle
HRX: ½ circle, round body needle - extra strong						
HRX 27 	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0430 DX 0431 DX 0432	
HRX 35 	3.5 4	0 1	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0281 DX 0282	
	3.5	0	1 x 90 cm	Violet	DX 0381	
HRX 37 	4	1	1 x 90 cm	Violet	DX 0384	

POLYDOX monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
 HRX 40	3.5	0	1 x 75 cm	Violet	DX 0285	
	3.5	0	1 x 90 cm	Violet	DX 0385	
	4	1	1 x 90 cm	Violet	DX 0386	
	5	2	1 x 90 cm	Violet	DX 0387	
	3.5	0	1 x 150 cm	Violet	DX 4035	Loop
	4	1	1 x 150 cm	Violet	DX 4040	Loop
	4	1	1 x 240 cm	Violet	DX 4024	Loop
 HRX 45	4	1	1 x 90 cm	Violet	DX 0390	
	5	2	1 x 90 cm	Violet	DX 0391	
 HRX 48	3.5	0	1 x 90 cm	Violet	DX 0400	
	3.5	0	1 x 150 cm	Violet	DX 4135	Loop
	4	1	1 x 150 cm	Violet	DX 4140	Loop
 HRX 50	3.5	0	1 x 90 cm	Violet	DX 0394	
	4	1	1 x 90 cm	Violet	DX 0395	
	5	2	1 x 90 cm	Violet	DX 0396	
	3.5	0	1 x 150 cm	Violet	DX 1392	Loop
	3.5	0	1 x 150 cm	Violet	DX 4535	Loop
	4	1	1 x 150 cm	Violet	DX 4540	Loop
	5	2	1 x 150 cm	Violet	DX 4550	Loop
	4	1	1 x 200 cm	Violet	DX 4542	Loop
 HRX 70	4	1	1 x 150 cm	Violet	DX 4543	Loop



POLYDOX monofilament

UK07
126

Absorbable - Polydiox monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 25 	2	3/0	1 x 75 cm	Violet	DX 0581	
HRT 26 	3	2/0	1 x 75 cm	Violet	DX 0572	
HRT 37 	4	1	1 x 90 cm	Violet	DX 0599	
HRT 48 	4	1	1 x 90 cm	Violet	DX 0608	
HRTX: ½ circle, round body needle – tapercut, extra strong						
HRTX 40 	4	1	1 x 90 cm	Violet	DX 0610	
HRF: ½ circle, flattened body needle, taper point						
HRF 26 	2 3 3.5	3/0 2/0 0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 9266 DX 9267 DX 9268	
HRF 31 	3	2/0	1 x 75 cm	Violet	DX 9271	

POLYDOX monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DR: 3/8 circle, round body needle						
2 x DR 10 	0.7	6/0	1 x 75 cm	Violet	DX 1120	
DR 12 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0120 DX 0121	
DR 17 	1.5	4/0	1 x 75 cm	Violet	DX 0132	
DR 20 	1.5 3	4/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0141 DX 0143	
DRT: 3/8 circle, round body needle - tapercut						
DRT 13 	1	5/0	1 x 75 cm	Violet	DX 0521	
FR: 5/8 circle, round body needle						
2 x FR 27 	3	2/0	1 x 75 cm	Violet	DX 5367	
FR 65 	3.5	0	1 x 150 cm	Violet	DX 4635	Loop
GR: Straight, round body needle						
2 x GR 60 	2	3/0	1 x 75 cm	Violet	DX 5396	
KR: Ski, round body needle						
KR 27 	2	3/0	1 x 40 cm	Violet	DX 4630	Loop
JRT: Hook shaped, round body needle - tapercut						
JRT 30 	4	1	1 x 75 cm	Violet	DX 0694	

POLYDOX monofilament

UK07
128

Absorbable - Polydiox monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
HS 18 	1.5	4/0	1 x 75 cm	Violet	DX 0155	
HS 21 	2	3/0	1 x 75 cm	Violet	DX 0166	
HS 22 	2 3	3/0 2/0	1 x 45 cm 1 x 45 cm	Violet Violet	DX 0167 DX 0170	
HS 27 	3	2/0	1 x 75 cm	Violet	DX 0168	
HS 30 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0171 DX 0172	
HS 37 	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0177 DX 0178 DX 0179	
	3.5	0	1 x 90 cm	Violet	DX 1178	
HS 40 	4	1	1 x 75 cm	Violet	DX 0180	
HS 48 	4	1	1 x 90 cm	Violet	DX 0184	

POLYDOX monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 DS 15	1	5/0	1 x 45 cm	Violet	DX 0190	
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0196 DX 0197	
 DS 16	1	5/0	1 x 45 cm	Violet	DX 0198	
 DS 19	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Violet Violet Violet	DX 0210 DX 0211 DX 0212	
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0215 DX 0216	
 DS 25	2 3	3/0 2/0	1 x 45 cm 1 x 45 cm	Violet Violet	DX 0222 DX 0223	
	1.5 2 3 3.5 4	4/0 3/0 2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet Violet Violet	DX 0226 DX 0227 DX 0228 DX 0229 DX 0230	
 DS 30	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0232 DX 0234 DX 0235	
DSC: 3/8 circle, conventional cutting needle						
 DSC 25	3	2/0	1 x 75 cm	Violet	DX 8228	




POLYDOX monofilament

UK07
130

Absorbable - Polydiox monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Violet Violet	DX 0201 DX 0202	Advanced cutting needle
DSM 15 	1.5	4/0	1 x 45 cm	Violet	DX 0206	Advanced cutting needle
DSM 19 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0207 DX 0208	Advanced cutting needle
GS: Straight, reverse cutting needle						
2 x GS 65 	2	3/0	1 x 75 cm	Violet	DX 3242	
Minimal invasive surgery						
	3.5 4	0 1	1 x 70 cm 1 x 70 cm	Violet Violet	MS 0134 MS 0135	Loop Loop



The CTI group

A constellation of 3 European manufacturers of high quality products for the healthcare sector

The CTI group consists of: CHIRANA T. Injecta in Slovakia, CHIRANA T. Injecta in the Czech Republic and Vitrex Medical in Denmark.

The intention of CTI group is to serve as an "umbrella-organization", where the individual companies collaborate to be "Committed to Innovation" and reach the joint goal:
- to develop, produce and distribute high quality, healthcare products.





UK07
219

The CTI group

CHIRANA T. Injecta, Slovakia

CHIRANA T. Injecta is a global, state of the art manufacturer of healthcare products, which was founded in 1935.

In 1947 CHIRANA T. Injecta started manufacturing glass syringes and they have ever since focused on development and production of high quality, disposable products for the healthcare sector. Through continuously updating production facilities and expanding their product range, they are today supplying:

- Hypodermic syringes
- Perfusion syringes
- Pen and Blood collection needles
- Insulin syringes
- Cannula and lancet needles

CHIRANA T. Injecta, Czech Republic

CHIRANA T. Injecta was founded in 1930. Since the early 1990s they have made substantial investments in upgrading the production resources, enabling the company to offer products in an unique quality.

Their more than 85 years of experience with manufacturing sutures has resulted in top priority within quality control, product development and sutures manufactured in line with state of the art technology. Their present range consist of:

- Surgical Sutures, Ligatures and Meshes
- Veterinary bottles and Eye needles
- Microsurgery and Cardiovascular sutures
- Antibacterial sutures

VITREX Medical, Denmark

Vitrex Medical was founded in 1965 and has since the 1970s mainly focused on supplying to the laboratory and analytical sector.

With more than 50 years of experience Vitrex Medical has become the largest manufacturer of capillary tubes in the world today by developing equipment, establishing technologies, and optimizing production processes. Vitrex Medical currently supplies:

- Capillary tubes and Micro-Hematocrit tubes
- Micro-pipettes
- Safety lancets
- Vacuum tubes and Blood collection needles
- Pen needles and Diabetic lancets



Version UK07

Designed by Vitrex Medical A/S

2023-02 Rev.007

CHIRANA T. Injecta, s.r.o.

**Komoranska 2148
143 00 Prague
Czech Republic**

www.cti-group.com