



Vitrex Alcohol Swabs

Vitrex Alcohol swabs are an easy and convenient way to cleanse the skin and help fight bacteria at the site of application. Use as an antiseptic skin cleanser before injections and blood sampling.

The size of the Disinfection Wipe is 30x60 mm. Each Alcohol Swab is individually sealed and disposable for convenience and safety.

Technical Specifications

Description	Size	Content	Pcs./Box	Article No.
Alcohol Swabs	30 x 60 mm	82% Ethanol	100	521213
Alcohol Swabs	30 x 60 mm	82% Ethanol + 0.5% Chlorhexidine	100	522213

Shelf life 5 years



Transit packaging Pallet

Packaging:
Primary packaging: 100 pcs.
Transit packaging: 10,000 pcs.
Pallet: 400,000 pcs.