



Vitrex Capillary & Melting Point Tubes

Vitrex Capillary Tubes are perfect for the handling and transport of small volumes of liquid, melting point determinations and many other purposes.

Melting Point tubes are uniform and cut at right angles. The tubes are packed as bulk in boxes. The tubes are also possible to pack in vials of 100 pcs. upon special request.

Technical Specifications

d (mm) +/- 0.05	D (mm) +/- 0.05	L (mm) +/- 1	Pcs.	Capillary Tubs Both ends open	Melting Point One end Closed
0.80	1.00	80	2000	210013	218013
0.80	1.00	100	2000	210113	218113
0.90	1.20	70	1000	210213	218213
0.90	1.20	100	1000	210313	218313
1.00	1.40	75	1000	210413	218413
1.00	1.40	80	1000	210513	218513
1.00	1.40	100	1000	210613	218613
1.00	1.40	125	1000	210713	218713
1.00	1.40	150	1000	210813	218813
1.00	1.40	200	5000	210913	218913
1.15	1.55	80	1000	211013	219013
1.15	1.55	100	1000	211113	219113
1.20	1.65	75	500	211213	219213
1.20	1.65	100	500	211313	219313
1.20	1.65	125	1000	211413	219413
1.20	1.65	150	1000	211513	219513
1.20	1.65	200	3000	211613	219613
1.40	1.90	125	500	211713	219713
1.50	2.00	75	500	211813	219813
1.50	2.00	100	500	211913	219913

For other dimensions, tolerances, sealing, impregnation or packaging, please contact us.

Shelf life 4 years

