



Vitrex End-to-End Pipettes

Vitrex End-to-End Pipettes are intended for applications where accurate capillary tube volume is required. The pipettes are in conformity with ISO 7550.

The pipettes are packed in glass vials with a dispenser top, that makes it easy to shake out one tube at a time and thereby ensuring hygienic handling.

The dispenser tops are colour-coded to indicate the tube volume in conformity with ISO 1769. Ten vials are packaged in a sturdy carton in order to protect the tubes from damage during handling or transport.


Vitrex tolerances

Vol. (μ l)	Accuracy	CV
0.5 - 10	$\pm 0.50\%$	0.75%
20 - 100	$\pm 0.50\%$	0.60%



Accessories

Product	Article No.
Pipette holder - Standard	140413
Pipette holder - Compact	140313





Technical Specifications

Colour Code	Vol. (µl)	L (mm) +/- 1mm	Pcs./box	Plain	Natrium Hep.	Ammon. Hep.	K2EDTA
○	0.5	15	10x100	177613	*	*	*
●	1	29	10x100	165013	167513 *	169013 *	*
●	2	29	10x100	165113	167613 *	169113 *	*
●	3	29	10x100	165213	167713 *	169213 *	*
●	4	29	10x100	165313	167813 *	169313 *	*
○	5	29	10x100	165413	167913 *	169413 *	170113 *
●	10	29	10x100	174313	177113 *	179213 *	174113 *
●	20	29	10x100	174513	177213 *	179313 *	174413 *
○	25	29	10x100	174613	177313 *	179413 *	178213 *
●	40	29	10x100	174813	177413 *	179513 *	*
●	44.7	29	10x50	175013	177713 *	178613 *	*
●	50	29	10x50	174913	177513 *	179613 *	171113 *

* Available upon request
For lithium hep. please contact us

Shelf life 4 years

