

NEW

Vitrex Safety Pen Needles With Dual Automatic Protective Shields

Vitrex Safety Pen Needles are intended for a safe and simple insulin injection via most insulin pens. The application is the same as normal insulin pens, but they offer an automatic double protection without extra operations after using. The coloured indicator shields in the front and back-end provide a visual, physical and audible confirmation of safety. The needles are sterile and individually packed.

SAFETY FUNCTION

Coloured automatic safety shields on front and back-end providing double protection after using.

COMFORTABLE INJECTION

Triple sharpened needle and top lubrication technology, which offers minimum penetration force.






SIMPLE AND SAFE APPLICATION

Easy to trigger the safety function after injection. Same application as normal pen needles, but with a safety function both blocking top and bottom needle after use.

VERSATILE


Uniform hub screw which fits most commonly used insulin pens.

Technical Specifications


Colour Code		Specification	Article No.
Purple		30G x 5 mm	531113
Dark blue		30G x 8 mm	531213
Yellow		31G x 5 mm	531313
Orange		31G x 6 mm	531413
Green		32G x 4 mm	531513

Shelf life 5 years

Before use:



After use:



Packaging:
Primary packaging: 100 pcs.
Transit packaging: 800 pcs.
Pallet: 57,600 pcs.