



Coming soon

Vitrex® Filter Volume Capillary Tubes

Vitrex Filter Volume Capillary Tubes are Mylar wrapped glass tubes with a hydrophobic filter preset to achieve precise volume delivery. The filter controls the volume and serves as a piston which dispenses the liquid when pressed downward by the plunger. The Mylar wrap provides enhanced safety with no compromise in accuracy, as it will help contain the liquid and protect against cuts from broken glass in the event of breakage. Vitrex Filter Volume Capillary Tubes offer a two-in-one solution as they can collect and dispense. The tubes are available plain and heparinised. They are packed in vials containing 50 tubes and 50 plungers or 100 tubes and 5 plungers. Other packings are available upon request.

Vitrex® Plastic Bulb Dispenser

Vitrex Plastic Bulb Dispenser is a simple, one-piece plastic tube and dispenser in one, that offers easy and accurate volume measurement. The dispenser is intended for collecting and dispensing micro samples of blood. The possibility of contamination with Vitrex Plastic Bulb Dispenser is eliminated due to breakage with no compromise in accuracy. The dispenser is available in volumes of 5 µl to 80 µl and packed in bags of 10 or 50 pcs.




Coming soon