



Chiravac™ Blood Collection needles

Chiravac Blood Collection needles are made of polypropylene hub and stainless steel tube. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The needles are intended for blood collection with vacuum tubes. They are individually packed in boxes of 100 pcs, cartons of 1000 pcs.

USER CONVENIENCE: Colour-coded needle hub, screw-on needle hub compatible with all available holders, easy-to-open packaging and thin wall tube material.

COMFORT: Grinded, precision lancet point and silicone treated tip for smooth penetration. Comfortable collection experience.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Colour Code	Size Metric	Size G ¹ / ₂	Article No.
Blue ●	0.60 x 38 mm	23G x 1 1/2"	CHBCN23112
Black ●	0.70 x 38 mm	22G x 1 1/2"	CHBCN22112
Green ●	0.80 x 38 mm	21G x 1 1/2"	CHBCN21112
Yellow ●	0.90 x 38 mm	20G x 1 1/2"	CHBCN20112
Pink ●	1.20 x 38 mm	18G x 1 1/2"	CHBCN18112
Blue ●	0.60 x 25 mm	23G x 1"	CHBCN23100
Black ●	0.70 x 25 mm	22G x 1"	CHBCN22100
Green ●	0.80 x 25 mm	21G x 1"	CHBCN21100
Yellow ●	0.90 x 25 mm	20G x 1"	CHBCN20100
Pink ●	1.20 x 25 mm	18G x 1"	CHBCN18100

Shelf life 5 years

