



Medoject® Safety Hypodermic needles

Medoject Safety Hypodermic needles have an integrated safety feature for one-handed manipulation.

| Colour Code | Size Metric | Size G/” | Article No. |
|-------------|--------------|--------------|-------------|
| Blue ● | 0.60 x 30 mm | 23G x 1 1/4” | CHS23114 |
| Black ● | 0.70 x 40 mm | 22G x 1 1/2” | CHS22112 |
| Green ● | 0.80 x 40 mm | 21G x 1 1/2” | CHS21112 |
| Yellow ● | 0.90 x 40 mm | 20G x 1 1/2” | CHS20112 |

Shelf life 5 years

Medoject® Blunt fill needles

Medoject Blunt fill needles are also available with filter of 5µm for aspiration filtration procedure, which risks particle elimination during procedure.

| Colour Code | Size Metric | Size G/” | Article No. |
|-------------|-----------------------------------|--------------|-------------|
| Red ● | 1.20 x 40 mm Blunt fill | 18G x 1 1/2” | CH18112BF |
| Pink ● | 1.20 x 40 mm Blunt fill w/ filter | 18G x 1 1/2” | CH18112F |
| Pink ● | 1.20 x 50 mm Blunt fill w/ filter | 18G x 2” | CH18200F |

Shelf life 5 years



According to EN, ISO standards.
Quality management certified according to EN ISO 13485:2016. CE marked.

Medoject® Special Hypodermic needles

Medoject Special Hypodermic needles are intended for a wide special application range. There is possibility of short bevel lancet point grinding.

| Colour Code | Size Metric | Size G/” | Article No. |
|--------------|---------------|--------------|-------------|
| Blue ● | 0.60 x 80 mm | 23G x 3 1/8” | CH23318 |
| Green ● | 0.80 x 100 mm | 21G x 4” | CH21400 |
| Green ● | 0.80 x 120 mm | 21G x 4 3/4” | CH21434 |
| Yellow ● | 0.90 x 70 mm | 20G x 2 4/5” | CH20245 |
| Pink ● | 1.20 x 70 mm | 18G x 2 4/5” | CH18245 |
| Pink ● | 1.20 x 100 mm | 18G x 4” | CH18400 |
| White ○ | 1.60 x 70 mm | 16G x 2 4/5” | CH16245 |
| Blue grey ● | 1.80 x 70 mm | 15G x 2 4/5” | CH15245 |
| Pale green ● | 2.00 x 70 mm | 14G x 2 4/5” | CH14245 |
| Pale green ● | 2.00 x 80 mm | 14G x 3 1/8” | CH14318 |
| Pale green ● | 2.00 x 100 mm | 14G x 4” | CH14400 |
| Pale green ● | 2.00 x 120 mm | 14G x 4 3/4” | CH14434 |

Other sizes available upon request.

Shelf life 5 years

Inox® Ophthalmic needles

Inox Ophthalmic needles has a blunt end for ophthalmic application, which includes scratch-free edges for comfort use..

| Colour Code | Size Metric | Size G/” | Article No. STRAIGHT | Article No. BENT |
|-------------|--------------|------------|----------------------|------------------|
| Grey ● | 0.40 x 20 mm | 27G x 4/5” | CHOP27045S | CHOP27045B |
| Blue ● | 0.60 x 20 mm | 23G x 4/5” | CHOP23045S | CHOP23045B |

Shelf life 5 years