



## Vitrex Sedimentation Systems

Vitrex Vacusence Sedimentation Tubes contain 3.2% sodium citrate for ESR test. The tubes are available with an open or closed sedimentation system.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Vacusence+ Article No.	Vacusence Article No.
ESR Tube	2	13 x 75	4NC/Sodium Citrate 3.2%	Paper	111602+	-

**Shelf life 14 months**

### Vacusence+



Transit packaging

Pallet

#### Packaging:

Primary packaging: 100 pcs.  
 Transit packaging: 1,200 pcs.  
 Pallet: 54,000 pcs.