



Chiraflex® IV Cannulas and Chiraflex® Safety IV Cannulas

Chiraflex IV Cannulas and Chiraflex Safety IV Cannulas are intended for peripheral vascular application. The catheter is made from Teflon material and available in Safety. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. They are packed in boxes of 50 pcs, cartons of 1000 pcs.

USER CONVENIENCE: Compatible with all standard medical devices for wide application.

USER SAFETY: Safety cannula is available.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

| Colour Code | Size G ¹ | Size Metric | Article No. With Port With Wings | Article No. Without Port With Wings | Article No. Without Port Without Wings | Article No. SAFETY With Port With Wings |
|-------------|---------------------|-------------|----------------------------------|-------------------------------------|--|---|
| Purple ● | 26G | 0.60 x 19 | CHFX0102619 | CHFX0202619 | CHFX0302619 | CHFXS0102619 |
| Yellow ● | 24G | 0.70 x 19 | CHFX0102419 | CHFX0202419 | CHFX0302419 | CHFXS0102419 |
| Blue ● | 22G | 0.90 x 25 | CHFX0102225 | CHFX0202225 | CHFX0302225 | CHFXS0102225 |
| Pink ● | 20G | 1.10 x 32 | CHFX0102032 | CHFX0202032 | CHFX0302032 | CHFXS0102032 |
| Green ● | 18G | 1.30 x 45 | CHFX0101845 | CHFX0201845 | CHFX0301845 | CHFXS0101845 |
| White ○ | 17G | 1.50 x 45 | CHFX0101745 | CHFX0201745 | CHFX0301745 | CHFXS0101745 |
| Grey ● | 16G | 1.70 x 45 | CHFX0101645 | CHFX0201645 | CHFX0301645 | CHFXS0101645 |
| Orange ● | 14G | 2.10 x 45 | CHFX0101445 | CHFX0201445 | CHFX0301445 | CHFXS0101445 |

Shelf life 5 years