



## Vitrex Soft Blood Glucose Lancets

Vitrex Soft Blood Glucose Lancets are intended for blood glucose measurement. The soft lancets are universal and compatible with most lancing devices especially Vitrex Compact Eject. The range includes standard soft lancets and siliconized soft lancets available in different sizes.

### GOOD BLOOD FLOW

Triple ground needle point.

### VERSATILE

Fit most commonly used lancing devices including the Vitrex Compact Eject lancing device.

### STERILE

Sterilised by gamma irradiation at minimum 25kGy.

### CONVENIENT

Boxes have dispensing feature.

### NEW SILICONIZED NEEDLE

For less pain and smoother penetration.

### Technical Specifications

Colour Code	Specification	Pcs./Box	Article No.
Light blue	23G (0.64 mm)	200	351213
Red	28G (0.36 mm)	200	352213
Green	30G (0.30 mm)	200	354213
Dark blue	32G (0.26 mm)	200	357213

**Shelf life 5 years**

### Technical Specifications Siliconized Needles

Colour Code	Specification	Pcs./Box	Article No.
Red	28G (0.36 mm)	200	352413
Green	30G (0.30 mm)	200	354413

**Shelf life 5 years**



Transit packaging

Pallet

#### Packaging:

Primary packaging: 200 pcs.  
Secondary packaging: 2,000 pcs.  
Transit packaging: 20,000 pcs.  
Pallet: 480,000 pcs.