

Silon

**Silon
monofilament**


Silon monofilament

Silon

monofila

Silon

monofilament



SILON monofilament

SILON monofilament is a non-absorbable, synthetic suture made of polyamide, which has a smooth non-porous surface.

The thread is extremely useful for skin sutures and has a high tensile strength in knots, and knots tie down easily. The thread is available in blue and black.

Silon monofilament can cause a slight initial inflammatory tissue reaction.



UK07
155

Indication: Plastic and General surgery and wherever suturing the skin.

Contraindication: Silon monofilament should not be used when long-term mechanical support is required.

Chemical Composition: Polyamide 6

Sizes: 10/0 to 2 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue and Black



SILON monofilament blue ligatures




UK07
156

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	405
1.5	4/0	3 x 45 cm	406
2	3/0	3 x 45 cm	407
3	2/0	3 x 45 cm	408
1	5/0	10 x 45 cm	409
1.5	4/0	10 x 45 cm	410
2	3/0	10 x 45 cm	411
3	2/0	10 x 45 cm	412
1	5/0	5 x 75 cm	413
1.5	4/0	5 x 75 cm	414












Pre-cut threads

EP	USP	Length	Cat. No
2	3/0	5 x 75 cm	415
3	2/0	5 x 75 cm	416
2	3/0	2 x 100 cm	1485
3.5	0	2 x 100 cm	1487
1	5/0	1 x 150 cm	401
1.5	4/0	1 x 150 cm	402
2	3/0	1 x 150 cm	403
3	2/0	1 x 150 cm	404
3.5	0	1 x 150 cm	1405
4	1	1 x 150 cm	1406



Non-absorbable - Silon monofilament










SILON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	1	5/0	1 x 45 cm	Blue	SM 2453	
HR 13 	1.5	4/0	1 x 45 cm	Blue	4125	
	0.7	6/0	1 x 60 cm	Blue	SM 4126	
HR 18 	1.5	4/0	1 x 75 cm	Blue	4104	
HR 22 	2	3/0	1 x 75 cm	Blue	4106	
HR 27 	3	2/0	1 x 75 cm	Blue	SM 4137	
HR 30 	3	2/0	1 x 75 cm	Blue	SM 4142	
HR 35 	3.5	0	1 x 75 cm	Blue	SM 4149	
HR 40 	4	1	1 x 75 cm	Blue	4152	
	3.5	0	1 x 100 cm	Blue	SM 4154	
HRX: ½ circle, round body needle - extra strong						
HRX 35 	3.5	0	1 x 75 cm	Blue	SM 4241	
HRX 40 	4	1	1 x 150 cm	Blue	SM 5033	Loop
HRX 45 	4	1	1 x 200 cm	Blue	SM 4540	Loop











SILON monofilament

UK07
158

Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 48 	4	1	1 x 150 cm	Blue	SM 5034	Loop
	4	1	1 x 200 cm	Blue	SM 5037	Loop
HRX 50 	3.5	0	1 x 150 cm	Blue	SM 5035	Loop
	4	1	1 x 200 cm	Blue	SM 5040	Loop
HRT: ½ circle, round body needle - tapercut						
HRT 13 	0.7 1 1.5	6/0 5/0 4/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Blue Blue Blue	SM 2552 SM 2553 SM 2554	
DR: 3/8 circle, round body needle						
DR 8 	0.5 1	7/0 5/0	1 x 45 cm 1 x 45 cm	Blue Blue	SM 2448 450	
DR 10 	1	5/0	1 x 45 cm	Blue	451	
DR 15 	1.5	4/0	1 x 45 cm	Blue	4024	
DR 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Blue Blue	4003 4021	
DR 21 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Blue Blue	SM 2325 SM 2326	
DR 31 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Blue Blue	SM 2330 SM 2331	

SILON monofilament







Needle	EP	USP	Length	Colour	Cat. No	Other information
GR: Straight, round body needle						
2 x GR 60 	3.5	0	1 x 90 cm	Blue	SM 2358	
HS: 1/2 circle, reverse cutting needle						
HS 12 	1	5/0	1 x 45 cm	Blue	SM 2305	
	1	5/0	1 x 75 cm	Blue	SM 2306	
HS 15 	2	3/0	1 x 45 cm	Blue	SM 2308	
HS 18 	1.5	4/0	1 x 75 cm	Blue	4420	
HS 25 	2	3/0	1 x 75 cm	Blue	SM 2311	
HS 40 	4	1	1 x 75 cm	Blue	SM 2413	
	3.5	0	1 x 100 cm	Blue	SM 3382	
DS: 3/8 circle, reverse cutting needle						
DS 10 	0.5 0.7	7/0 6/0	1 x 45 cm 1 x 45 cm	Blue Blue	SM 2251 SM 2252	
DS 11 	0.5	7/0	1 x 45 cm	Blue	SM 2254	
DS 12 	0.5 0.7 0.7 1 1 1.5	7/0 6/0 6/0 5/0 5/0 4/0	1 x 45 cm 1 x 45 cm 1 x 45 cm 1 x 45 cm 1 x 45 cm 1 x 45 cm	Blue Blue Black Blue Black Blue	4328 SM 2255 SM 4255 SM 2256 SM 4256 SM 2257	
DS 13 	1	5/0	1 x 75 cm	Blue	SM 2248	










SILON monofilament

UK07
160

Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 DS 15	0.7	6/0	1 x 45 cm	Blue	SM 2260	
	1	5/0	1 x 45 cm	Blue	SM 2261	
	1	5/0	1 x 45 cm	Black	SM 4261	
	1.5	4/0	1 x 45 cm	Blue	SM 2262	
	2	3/0	1 x 45 cm	Blue	SM 2263	
	1.5	4/0	1 x 75 cm	Blue	SM 2362	
 DS 16	0.7	6/0	1 x 45 cm	Blue	SM 2460	
	0.7	6/0	1 x 60 cm	Blue	SM 2459	
 DS 19	0.7	6/0	1 x 45 cm	Blue	4332	
	1	5/0	1 x 45 cm	Blue	SM 2266	
	1	5/0	1 x 45 cm	Black	SM 4266	
	1.5	4/0	1 x 45 cm	Blue	SM 2367	
	2	3/0	1 x 45 cm	Blue	SM 2368	
	2	3/0	1 x 45 cm	Black	SM 4368	
	3	2/0	1 x 45 cm	Blue	SM 2369	
	1	5/0	1 x 75 cm	Blue	SM 2366	
	1.5	4/0	1 x 75 cm	Blue	SM 2267	
	2	3/0	1 x 75 cm	Blue	SM 2268	
3	2/0	1 x 75 cm	Blue	SM 2269		
	1	5/0	1 x 90 cm	Blue	SM 2466	
 DS 20	1.5	4/0	1 x 75 cm	Blue	SM 2471	
 DS 24	1.5	4/0	1 x 75 cm	Blue	SM 3277	
	2	3/0	1 x 75 cm	Blue	SM 3276	
	3	2/0	1 x 75 cm	Blue	SM 3275	
 DS 25	1.5	4/0	1 x 45 cm	Blue	SM 2372	
	1.5	4/0	1 x 45 cm	Black	SM 4372	
	2	3/0	1 x 45 cm	Blue	SM 2373	
	2	3/0	1 x 45 cm	Black	SM 4373	
	3	2/0	1 x 45 cm	Blue	SM 2374	
	3	2/0	1 x 45 cm	Black	SM 4374	
	1	5/0	1 x 75 cm	Blue	SM 2375	
	1.5	4/0	1 x 75 cm	Blue	SM 2273	
	2	3/0	1 x 75 cm	Blue	SM 2270	
	3	2/0	1 x 75 cm	Blue	SM 2271	
3.5	0	1 x 75 cm	Blue	SM 2272		
4	1	1 x 75 cm	Blue	SM 2274		
	2	3/0	1 x 90 cm	Blue	SM 3272	
	3	2/0	1 x 90 cm	Blue	SM 3271	

SILON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 <p>DS 30</p>	3	2/0	1 x 45 cm	Blue	SM 2384	
	3.5	0	1 x 45 cm	Blue	SM 2385	
	1.5	4/0	1 x 75 cm	Blue	SM 4324	
	2	3/0	1 x 75 cm	Blue	4325	
 <p>DS 37</p>	3	2/0	1 x 75 cm	Blue	SM 2284	
	3.5	0	1 x 75 cm	Blue	SM 2285	
	2	3/0	1 x 90 cm	Blue	SM 4325	
	3	2/0	1 x 90 cm	Blue	SM 4326	
 <p>DS 38</p>	4	1	1 x 90 cm	Blue	SM 4327	
	3.5	0	1 x 75 cm	Blue	SM 2290	
	2	3/0	1 x 75 cm	Blue	SM 2291	
	3	2/0	1 x 75 cm	Blue	SM 2292	
 <p>DS 40</p>	4	1	1 x 75 cm	Blue	SM 2293	
	3	2/0	1 x 90 cm	Blue	SM 3291	
	3	2/0	1 x 90 cm	Blue	SM 2497	
	3.5	0	1 x 90 cm	Blue	SM 2495	
 <p>DS 50</p>	4	1	1 x 100 cm	Blue	SM 2693	
	4	1	1 x 100 cm	Blue	SM 2693	
 <p>DS 60</p>	3	2/0	1 x 100 cm	Blue	SM 3388	
	4	1	1 x 100 cm	Blue	SM 3390	
	4	1	1 x 100 cm	Blue	SM 3395	
	4	1	1 x 100 cm	Blue	SM 3395	
 <p>DS 90</p>	4	1	1 x 100 cm	Blue	SM 3395	
	4	1	1 x 100 cm	Blue	SM 3395	



SILON monofilament

UK07
162

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSC: 3/8 circle, conventional cutting needle						

DSC 90



3	2/0	1 x 100 cm	Blue	SM 8397
3.5	0	1 x 100 cm	Blue	SM 8394

DSM: 3/8 circle, reverse cutting needle - micropoint

DSM 12



0.7	6/0	1 x 45 cm	Blue	SM 2055
1	5/0	1 x 45 cm	Blue	SM 2056
1.5	4/0	1 x 45 cm	Blue	SM 2057
1.5	4/0	1 x 45 cm	Black	SM 4057
1	5/0	1 x 75 cm	Blue	SM 2053

Advanced cutting needle

DSM 15



1	5/0	1 x 45 cm	Blue	SM 2060
1.5	4/0	1 x 45 cm	Blue	SM 2061
2	3/0	1 x 45 cm	Blue	SM 2062

Advanced cutting needle

DSM 19



1	5/0	1 x 45 cm	Blue	SM 2064
1.5	4/0	1 x 45 cm	Blue	SM 2065
2	3/0	1 x 45 cm	Blue	SM 2066
1	5/0	1 x 75 cm	Blue	SM 2166
2	3/0	1 x 75 cm	Blue	SM 2168

Advanced cutting needle

DSM 25



1.5	4/0	1 x 45 cm	Blue	SM 2073
2	3/0	1 x 45 cm	Blue	SM 2074
2	3/0	1 x 75 cm	Blue	SM 2078
3	2/0	1 x 75 cm	Blue	SM 2079

Advanced cutting needle

DSM 40



3	2/0	1 x 75 cm	Blue	SM 2090
---	-----	-----------	------	---------

Advanced cutting needle

DSMC: 3/8 circle, conventional cutting needle - micropoint

DSMC 19



1	5/0	1 x 100 cm	Blue	SM 8081
---	-----	------------	------	---------

Advanced cutting needle




DSMC 20



1.5	4/0	1 x 45 cm	Blue	SM 8065
-----	-----	-----------	------	---------

Advanced cutting needle

SILON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
GS: Straight, reverse cutting needle						
GS 50 	3	2/0	1 x 90 cm	Blue	SM 3295	
GS 60 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Blue Blue	SM 3283 SM 3284	
GS 75 	3	2/0	1 x 90 cm	Blue	SM 3294	
	3.5	0	1 x 75 cm	Blue	SM 3299	
Minimal invasive surgery						
	3.5	0	1 x 70 cm	Blue	MS 1124	Loop

