



## Vitrex Blood Gas Tubes

Vitrex capillary tubes for blood gas analyses are made from heavy walled glass, suitable for 2.3 mm Ø adaptors. The uniform heparinisation secure proper function of the tubes for blood gas analysis. Vitrex Blood Gas Tubes are also available in plastic with sodium heparin in the same sizes as Vitrex Blood Gas Tubes in glass.

### Accessories

Stirrers packed in vials of 250 pcs.  
End Caps packed in boxes of 500 pcs.  
Magnets packed in bags of 10 pcs.





### Blood Gas Tubes

Vol. ( µl )	Pcs.	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 1 mm	Na Hep. Article No.	Li Hep. Article No.
130 (min)	4 x 250	1.60	2.30	75	182513	182613
175 (min)	4 x 250	1.60	2.30	100	182413	183313
220 (min)	4 x 250	1.60	2.30	125	182313	182013

Special types upon request.

**Shelf life 4 years**

### Plastic Blood Gas Tubes

Vol. ( µl )	Pcs.	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 1 mm	Plastic Na Hep. Article No.
130 (min)	4 x 200	1.60	2.30	75	282513
175 (min)	4 x 200	1.60	2.30	100	282413
220 (min)	4 x 200	1.60	2.30	125	282313

Other types available upon request.

**Shelf life 2 years**

### Accessories

Description	For tubes with D	Pcs.	Article No.
End caps for D (mm)	2.20 - 2.80 mm	10 x 500	140813
End caps for D (mm)	1.50 - 2.20 mm	10 x 500	140713
Stirrers		10 x 250	140913
Magnets		1 x 10	142400

