


Vitrex Precalibrated Pipettes

Vitrex Precalibrated Pipettes are intended for applications where accurate capillary tube volume is required. The pipettes are in conformity with ISO 7550.















The pipettes are packed in cardboard vials. The pipettes have a standard length of 125 mm and the colour-coding on each pipette indicates the volume, which is measured from the base and up to the thin black ringmark.

Vitrex tolerances

Vol. (μl)	Accuracy	CV
1-5	$\pm 0.30\%$	0.60%
10-200	$\pm 0.25\%$	0.50%



Technical Specifications

														
Vol. (µl)	1,2,3,4,5	5	10	5-10	20	25	40	20-40	44.7	50	20-50	100	50-100	100-200
Pcs.	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x125
Plain	128813	126013	126213	125513	126413	126613	129313	125913	126813	127013	125613	127213	125713	129413
Na Hep.	128913	127413	127613	125213	127813	128013	128413	129913	128213	128513	129113	128613	127113	127713

250 pcs. except 100-200 µl, which are packed in vials of 125 pcs.
Other coatings available upon request i.e. Lithium, Ammonium, EDTA.

Shelf life 4 years

Vol.