



Chirahem™ Transfusion set

Chirahem Transfusion set is intended for blood transfusion administration. The set has a 200 µm blood filter.

Product	Chamber Vol. (ml)	Tube length (mm)	Packaging Type	Pcs./Carton	Article No.
Chirahem	17	1500	Blister	100	CHT314
Chirahem Dehp free	17	1500	Blister	100	CHT514
Chirahem	17	1500	Poly bag	300	CHT320
Chirahem with needle	17	1500	Poly bag	300	CHT321
Chirahem with metal spike	17	1500	Poly bag	300	CHT322
Chirahem latex free	17	1500	Blister	100	CHT518

Shelf life 5 years